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Indications for Use: KeyPrint® KeySplint Soft™ is indicated for the fabrication of orthodontic and dental appliances such as bite planes, mouthguards, nightguards, snoring appliances, splints and repositionalers.

Product Description: KeyPrint® KeySplint Soft™ is a light-curing resin for the 3D printing of flexible biocompatible dental devices for use in DLP 3D printers utilizing wavelengths between 385nm – 405nm.

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Contraindications:

Contains methacrylate

monomers and oligomers which, although rare, may cause an allergic reaction in individuals sensitive to acrylic containing products.

Warnings & Precautions:

- Deviation from the described manufacturing process may compromise biocompatibility, user safety, and lead to unwanted material properties of the finished printed product.

2. Follow all recommended validated settings for your printer and post curebox located on the Keystone website to achieve biocompatible and effective print results.

3. Review the product Safety Data Sheet (SDS) prior to use.

4. As per the SDS, wear proper personal protective equipment when handling KeyPrint® resins and uncured printed parts.

5. Avoid contact with skin and eyes when handling this resin and uncured printed parts. In case of accidental contact, follow the "First-aid measures" listed in Section 4 of the SDS and seek medical attention if necessary.

6. Perform final processing (i.e., polishing).

9. Prior to delivery to the patient, clean the dental oral appliance with soap and water to ensure the device is free of any debris from the polishing process.

10. Part is ready for use. The finished medical device resulting from these instructions for use is safe, biocompatible and effective.

CAUTION:

Federal law restricts this device to sale by, or on the order of a dental professional.

Directions for Use:

- Ensure that resin is tempered to ambient temperature (20-25°C/68-77°F) prior to printing.

Note: While handling KeyPrint® KeySplint Soft™, we recommend wearing personal protective equipment (i.e., safety glasses, lab coat, closed-toe shoes, gloves, etc.).

2. In order to achieve consistency of the resin and to prevent bubbles, agitate the bottle 1 hour prior to use. If bubbles are present, remove with a clean instrument/spatula.

3. Only use KeyPrint® KeySplint Soft™ predetermined settings for your DLP 3D printer and post curebox provided on Keystone's website. KeySplint Soft™ should be used with a 385nm - 405nm UV light source. Printers using alternative light sources require validation by Keystone's technical team for optimal settings.

4. Resin coated printed parts should be cleaned of any residual liquid resin with Isopropyl Alcohol (IPA) (at least 97% within approximately 8 hours from the completion of the print).

Note: Do not allow the parts to sit in IPA for longer than 5 minutes, as the properties may begin to deteriorate. Keystone discourages the use of denatured alcohol or ethanol for cleaning as it may diminish or degrade the quality of the finished parts.

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