



Liqcreate Dental Model Pro Beige

Liqcreate Dental Model Pro Beige is perfect for manufacturing highly accurate dental models, dies and aligner models. This high accuracy beige colored resin is optimized for Mask Stereolithography (MSLA), Liquid Crystal Display (LCD), Digital Light Processing (DLP) and Stereolithography (SLA) 3D-printers. This resin has been developed and validated in collaboration with dental companies and professionals in order to meet the high expectations in this disruptive market.

Product description

Liqcreate Dental Model Pro Beige is a matte, opaque beige colored photopolymer. Its matte surface finish accentuates depth and detail in dental models, which accommodates a perfect view on undercuts and other important markers. In addition, 3D-printed models from Dental Model Pro Beige are excellent for crown & bridge build-ups. 3D-printed parts from this material have exceptional dimensional stability and extremely low shrinkage during printing. Liqcreate Dental Model Pro Beige is easy to use on all open source SLA, MSLA and DLP 3D-printers in the range of 385 - 420nm. This material has excellent properties like low shrinkage and low odor, accuracy and dimensional stability, making it perfect for the production of dental models, C&B, thermoforming and aligner models.

Key benefits

- Excellent dimensional stability
- High detail & matte surface
- Low shrinkage
- Low Odor


3D-Printer compatibility

- Asiga Max & 4K
- Ackuretta Dentiq
- Photon and Phrozen Series
- Open 385 - 420nm LCD/MSLA, DLP and SLA 3D-printers

Order information

Order directly at the [Liqcreate store](https://www.liqcreate.com) or send your inquiry to order@liqcreate.com with the following order numbers.

Liqcreate Dental Model Pro Beige	250gram	Order number LDMP00250
Liqcreate Dental Model Pro Beige	1 kg	Order number LDMP01000





Liqcreate Dental Model Pro Beige Technical Data

Liquid properties			
Appearance	Opaque beige liquid	E _c	4.29 mJ/cm ²
Viscosity	700 cps at 25 °C	D _p metric	0.12 mm
Density	1.18 g/cm ³	D _p imperial	4.89 mils

Polymer properties			
Mechanical properties		30 minutes high power LED curing at 40 °C	
Description	ASTM Method	Metric	Imperial
Tensile Strength	D638M	37 MPa	5.4 ksi
Tensile Modulus	D638M	1.9 GPa	276 ksi
Elongation at break	D638M	4.8 %	
Flexural Strength	D790M	64 MPa	9.3 ksi
Flexural modulus	D2240	1.7 GPa	247 ksi
Shore D Hardness	D2240	85	
Water sorption	D570-98	0.37 %	
Linear shrinkage during printing	Internal method	0.3 %	
Linear shrinkage during UV-curing	Internal method	0.4 %	

These values may vary and depend on individual machine processing and post-curing. Liqcreate Dental Model Pro Beige is not tested on biocompatibility, therefore the material should not be used for applications inside the human body.

Visit www.liqcreate.com for more information about this product.